

For Research Use Only

ISG15 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68408-2-PBS



Basic Information

Catalog Number: 68408-2-PBS	GenBank Accession Number: BC009507	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 9636	CloneNo.: 3F1B1
Isotype: IgG1	UNIPROT ID: P05161	
Immunogen Catalog Number: AG8953	Full Name: ISG15 ubiquitin-like modifier	
	Calculated MW: 165 aa, 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51686-1, ISG15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68408-2-PBS. Detection antibody: 68408-3-PBS.



Cytometric bead array standard curve of MP51686-1, 1SG15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68408-2-PBS. Detection antibody: 68408-3-PBS. Standard: Ag8953. Range: 3.125-400 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 68408-2-Ig against Human ISG15 were performed. The affinity constant is 209 nM.