

For Research Use Only

ISG15 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68408-3-PBS



Basic Information

Catalog Number:

68408-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8953

GenBank Accession Number:

BC009507

GeneID (NCBI):

9636

UNIPROT ID:

P05161

Full Name:

ISG15 ubiquitin-like modifier

Calculated MW:

165 aa, 18 kDa

Purification Method:

Protein G purification

CloneNo.:

2A9C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51686-1, ISG15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68408-2-PBS. Detection antibody: 68408-3-PBS.



Cytometric bead array standard curve of MP51686-1, ISG15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68408-2-PBS. Detection antibody: 68408-3-PBS. Standard:Ag8953. Range: 3.125-400 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 68408-3-Ig against Human ISG15 were performed. The affinity constant is 734 nM.