

For Research Use Only

SYAP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68410-3-PBS



Basic Information

Catalog Number:

68410-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9626

GenBank Accession Number:

BC014657

GeneID (NCBI):

94056

UNIPROT ID:

Q96A49

Full Name:

synapse associated protein 1, SAP47 homolog (Drosophila)

Calculated MW:

352 aa, 40 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H4E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51310-1, SYAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68410-2-PBS. Detection antibody: 68410-3-PBS. Standard:Ag9626. Range: 0.098-100 ng/mL.