

For Research Use Only

# PNMA1 Monoclonal antibody, PBS Only

Catalog Number: 68415-1-PBS



## Basic Information

Catalog Number:

68415-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33143

GenBank Accession Number:

BC039577

GeneID (NCBI):

9240

UNIPROT ID:

Q8ND90

Full Name:

paraneoplastic antigen MA1

Calculated MW:

353 aa, 40 kDa

Observed MW:

37 kDa

Purification Method:

Protein G purification

CloneNo.:

2D8C6

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 68415-1-Ig (PNMA1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68415-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using PNMA1 antibody (68415-1-Ig, Clone: 2D8C6 ) at dilution of 1:800 and CoraLite®488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 68415-1-PBS in a different storage buffer formulation.