

For Research Use Only

# SSSCA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68423-2-PBS



## Basic Information

Catalog Number:

68423-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16783

GenBank Accession Number:

BC014791

GeneID (NCBI):

10534

UNIPROT ID:

O60232

Full Name:

Sjogren syndrome/scleroderma  
autoantigen 1

Calculated MW:

199 aa, 21 kDa

Purification Method:

Protein G purification

CloneNo.:

1D8H12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Cytometric bead array sample test of MP51566-1, SSSCA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68423-2-PBS. Detection antibody: 68423-3-PBS.



Cytometric bead array standard curve of MP51566-1, SSSCA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68423-2-PBS. Detection antibody: 68423-3-PBS. Standard:Ag16783. Range: 0.781-100 ng/mL.