

For Research Use Only

# ZNF189 Monoclonal antibody, PBS Only



Catalog Number: 68424-1-PBS

## Basic Information

Catalog Number:

68424-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG33164

GenBank Accession Number:

BC126120

GeneID (NCBI):

7743

UNIPROT ID:

O75820

Full Name:

zinc finger protein 189

Observed MW:

70-73 kDa

Purification Method:

Protein A purification

CloneNo.:

2A4B4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68424-1-Ig (ZNF189 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68424-1-PBS in a different storage buffer formulation.



HeLa cells were subjected to SDS PAGE followed by western blot with 68424-1-Ig (ZNF189 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68424-1-PBS in a different storage buffer formulation.