

For Research Use Only

ATP6V1H Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68425-2-PBS



Basic Information

Catalog Number:	68425-2-PBS	GenBank Accession Number:	BC025275	Purification Method:	Protein G purification
Source:	Mouse	GeneID (NCBI):	51606	CloneNo.:	2E6D7
Isotype:	IgG1	UNIPROT ID:	Q9UI12		
Immunogen Catalog Number:	AG10130	Full Name:	ATPase, H ⁺ transporting, lysosomal 50/57kDa, V1 subunit H		
		Calculated MW:	483 aa, 56 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51766-1, ATP6V1H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68425-2-PBS. Detection antibody: 68425-3-PBS.



Cytometric bead array standard curve of MP51766-1, ATP6V1H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68425-2-PBS. Detection antibody: 68425-3-PBS. Standard: Ag10130. Range: 0.781-100 ng/mL