

For Research Use Only

ATP6V1H Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68425-3-PBS



Basic Information

Catalog Number:

68425-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10130

GenBank Accession Number:

BC025275

GeneID (NCBI):

51606

UNIPROT ID:

Q9UI12

Full Name:

ATPase, H⁺ transporting, lysosomal
50/57kDa, V1 subunit H

Calculated MW:

483 aa, 56 kDa

Purification Method:

Protein G purification

CloneNo.:

2A4C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51766-1, ATP6V1H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68425-2-PBS. Detection antibody: 68425-3-PBS.



Cytometric bead array standard curve of MP51766-1, ATP6V1H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68425-2-PBS. Detection antibody: 68425-3-PBS. Standard: Ag10130. Range: 0.781-100 ng/mL