

For Research Use Only

BIK Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68438-4-PBS



Basic Information

Catalog Number:

68438-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6353

GenBank Accession Number:

BC001599

GeneID (NCBI):

638

UNIPROT ID:

Q13323

Full Name:

BCL2-interacting killer (apoptosis-inducing)

Calculated MW:

18 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D11D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50659-2, BIK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68438-4-PBS. Detection antibody: 68438-3-PBS. Standard:Ag6353. Range: 0.098-6.25 ng/mL



Cytometric bead array standard curve of MP50659-3, BIK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68438-4-PBS. Detection antibody: 68438-2-PBS. Standard:Ag6353. Range: 0.098-12.5 ng/mL