

For Research Use Only

CUX1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68449-4-PBS



Basic Information

Catalog Number:

68449-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28528

GenBank Accession Number:

BC025422

GeneID (NCBI):

1523

UNIPROT ID:

Q13948

Full Name:

cut-like homeobox 1

Calculated MW:

676 aa, 77 kDa

Purification Method:

Protein G purification

CloneNo.:

1F11C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51158-2, CUX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68449-4-PBS. Detection antibody: 68449-5-PBS. Standard:Ag28528. Range: 3.125-100 ng/mL