

For Research Use Only

FBLN1 Monoclonal antibody

Catalog Number: 68453-1-Ig



Basic Information

Catalog Number:

68453-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14249

GenBank Accession Number:

BC022497

GeneID (NCBI):

2192

UNIPROT ID:

P23142

Full Name:

fibulin 1

Calculated MW:

703 aa, 77 kDa

Observed MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

1F4A1

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : human placenta tissue, human plasma

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 68453-1-Ig (FBLN1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 68453-1-Ig (FBLN1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.