

For Research Use Only

SQRDL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68459-5-PBS



Basic Information

Catalog Number:

68459-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG10762

GenBank Accession Number:

BC016836

GeneID (NCBI):

58472

UNIPROT ID:

Q9Y6N5

Full Name:

sulfide quinone reductase-like (yeast)

Calculated MW:

450 aa, 50 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H7E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50545-3, SQRDL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68459-2-PBS. Detection antibody: 68459-5-PBS. Standard:Ag10762. Range: 1.563-100 ng/mL.