

For Research Use Only

# USP9X Monoclonal antibody, PBS Only

Catalog Number: 68460-1-PBS

Featured Product



## Basic Information

Catalog Number:

68460-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_001039590

GeneID (NCBI):

8239

UNIPROT ID:

Q93008

Full Name:

ubiquitin specific peptidase 9, X-linked

Calculated MW:

292 kDa

Observed MW:

260 kDa

Purification Method:

Protein G purification

CloneNo.:

1H6G11

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, pig

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68460-1-Ig (USP9X antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68460-1-PBS in a different storage buffer formulation.



WB result of USP9X antibody (68460-1-Ig; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-USP9X transfected HEK-293 cells. This data was developed using the same antibody clone with 68460-1-PBS in a different storage buffer formulation.