

For Research Use Only

LMNB2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68475-3-PBS



Basic Information

Catalog Number:

68475-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33296

GenBank Accession Number:

BC006551

GeneID (NCBI):

84823

UNIPROT ID:

Q03252

Full Name:

lamin B2

Calculated MW:

68 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E9A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sample test of MP51397-1, LMNB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68475-2-PBS. Detection antibody: 68475-3-PBS.



Cytometric bead array standard curve of MP51397-1, LMNB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68475-2-PBS. Detection antibody: 68475-3-PBS. Standard: Ag33296. Range: 0.195-100 ng/mL



Bi-layer interferometry (BLI) kinetic assays of 68475-3-Ig against Human LMNB2 were performed. The affinity constant is 50.7 nM.