

For Research Use Only

FLAD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68491-3-PBS



Basic Information

Catalog Number:

68491-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6845

GenBank Accession Number:

BC011378

GeneID (NCBI):

80308

UNIPROT ID:

Q8NFF5

Full Name:

FAD1 flavin adenine dinucleotide synthetase homolog (S. cerevisiae)

Calculated MW:

446 aa, 49 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D6C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50710-1, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-2-PBS. Detection antibody: 68491-3-PBS. Standard:Ag6845. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50710-2, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-4-PBS. Detection antibody: 68491-3-PBS. Standard:Ag6845. Range: 0.098-100 ng/mL