

For Research Use Only

# ATP1A2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68504-1-PBS



## Basic Information

Catalog Number:

68504-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG10489

GenBank Accession Number:

BC052271

GeneID (NCBI):

477

UNIPROT ID:

P50993

Full Name:

ATPase, Na<sup>+</sup>/K<sup>+</sup> transporting, alpha 2  
(+) polypeptide

Calculated MW:

1020 aa, 112 kDa

Purification Method:

Protein A purification

CloneNo.:

3D1C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50523-2, ATP1A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68505-2-PBS. Detection antibody: 68504-1-PBS. Standard:Ag10489. Range: 0.781-50 ng/mL