

For Research Use Only

TEX264 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68507-3-PBS



Basic Information

Catalog Number:

68507-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33710

GenBank Accession Number:

BC008742

GeneID (NCBI):

51368

UNIPROT ID:

Q9Y6I9

Full Name:

testis expressed 264

Calculated MW:

313 aa, 34 kDa

Purification Method:

Protein G purification

CloneNo.:

1D8H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50870-2, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-2-PBS. Detection antibody: 68507-3-PBS. Standard:Ag33710. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50870-3, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-4-PBS. Detection antibody: 68507-3-PBS. Standard:Ag33710. Range: 3.125-100 ng/mL.