

For Research Use Only

POLR2F Monoclonal antibody

Catalog Number: 68511-1-Ig



Basic Information

Catalog Number:

68511-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7472

GenBank Accession Number:

BC003582

GeneID (NCBI):

5435

UNIPROT ID:

P61218

Full Name:

polymerase (RNA) II (DNA directed)

polypeptide F

Calculated MW:

14 kDa

Observed MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

2D4D12

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : LNCaP cells, NIH/3T3 cells, MCF-7 cells, A549 cells, HeLa cells, 4T1 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68511-1-Ig (POLR2F antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 68511-1-Ig (POLR2F antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.