

For Research Use Only

GATM Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68513-4-PBS



Basic Information

Catalog Number: 68513-4-PBS	GenBank Accession Number: BC004141	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2628	CloneNo.: 3D5D10
Source: Mouse	UNIPROT ID: P50440	
Isotype: IgG1	Full Name: glycine amidinotransferase (L-arginine:glycine amidinotransferase)	
Immunogen Catalog Number: AG3754	Calculated MW: 423 aa, 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50565-2, GATM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68513-4-PBS. Detection antibody: 68513-5-PBS. Standard:Ag3754. Range: 1.563-25 ng/mL



Cytometric bead array standard curve of MP50565-3, GATM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68513-4-PBS. Detection antibody: 68513-6-PBS. Standard:Ag3754. Range: 0.781-12.5 ng/mL