

For Research Use Only

Cytokeratin 12 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68545-3-PBS



Basic Information

Catalog Number:

68545-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20332

GenBank Accession Number:

BC156641

GeneID (NCBI):

3859

UNIPROT ID:

Q99456

Full Name:

keratin 12

Calculated MW:

494 aa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E9A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51143-1, KRT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68545-3-PBS. Detection antibody: 68545-4-PBS. Standard:Ag20332. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP51143-3, KRT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68545-3-PBS. Detection antibody: 68545-7-PBS. Standard:Ag20332. Range: 0.098-1.563 ng/mL