

For Research Use Only

SUCLA2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68549-2-PBS

Basic Information

Catalog Number: 68549-2-PBS	GenBank Accession Number: BC027587	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 8803	CloneNo.: 2C11G2
Source: Mouse	UNIPROT ID: Q9P2R7	
Isotype: IgG1	Full Name: succinate-CoA ligase, ADP-forming, beta subunit	
Immunogen Catalog Number: AG17967	Calculated MW: 463 aa, 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50128-1, SUCLA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68549-2-PBS. Detection antibody: 68549-3-PBS. Standard:Ag17967. Range: 3.125-100 ng/mL.