

For Research Use Only

SUCLA2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68549-3-PBS

Basic Information

Catalog Number:

68549-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17967

GenBank Accession Number:

BC027587

GeneID (NCBI):

8803

UNIPROT ID:

Q9P2R7

Full Name:

succinate-CoA ligase, ADP-forming, beta subunit

Calculated MW:

463 aa, 50 kDa

Purification Method:

Protein G purification

CloneNo.:

3E5C2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50128-1, SUCLA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68549-2-PBS. Detection antibody: 68549-3-PBS. Standard:Ag17967. Range: 3.125-100 ng/mL.