

For Research Use Only

# TRPS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68562-3-PBS



## Basic Information

Catalog Number: 68562-3-PBS	GenBank Accession Number: NM_014112	Purification Method: Protein G purification
Concentration: 1000 µg/ml	GeneID (NCBI): 7227	CloneNo.: 1A10A1
Source: Mouse	UNIPROT ID: Q9UHF7	
Isotype: IgG1	Full Name: trichorhinophalangeal syndrome I	
Immunogen Catalog Number: AG30199	Calculated MW: 142 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50280-1, TRPS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68562-2-PBS. Detection antibody: 68562-3-PBS. Standard:Ag30199. Range: 0.781-200 ng/mL.