

For Research Use Only

TBC1D24 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68566-2-PBS



Basic Information

Catalog Number:

68566-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19380

GenBank Accession Number:

BC127015

GeneID (NCBI):

57465

UNIPROT ID:

Q9ULP9

Full Name:

TBC1 domain family, member 24

Calculated MW:

559 aa, 63 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E7A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51215-1, TBC1D24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68566-2-PBS. Detection antibody: 68566-3-PBS. Standard:Ag19380. Range: 3.125-200 ng/mL.