

For Research Use Only

TRADD Monoclonal antibody

Catalog Number: 68573-1-Ig



Basic Information

Catalog Number:

68573-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7768

GenBank Accession Number:

BC004491

GeneID (NCBI):

8717

UNIPROT ID:

Q15628

Full Name:

TNFRSF1A-associated via death domain

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Protein G purification

CloneNo.:

1C3H3

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, pig, rat

Positive Controls:

WB : A549 cells, LNCaP cells, HeLa cells, HEK-293 cells, Jurkat cells, K-562 cells, HSC-T6 cells, RAW264.7 cells, pig spleen tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68573-1-Ig (TRADD antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.