

For Research Use Only

SCNN1A Monoclonal antibody, PBS Only



Catalog Number: 68581-1-PBS

Basic Information

Catalog Number:

68581-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29548

GenBank Accession Number:

BC006526

GeneID (NCBI):

6337

UNIPROT ID:

P37088

Full Name:

sodium channel, nonvoltage-gated 1
alpha

Calculated MW:

76 kDa

Observed MW:

~80 kDa

Purification Method:

Protein A purification

CloneNo.:

1C2C2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68581-1-Ig (SCNN1A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68581-1-PBS in a different storage buffer formulation.