

For Research Use Only

TMEM132A Monoclonal antibody

Catalog Number: 68589-1-Ig



Basic Information

Catalog Number:

68589-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18149

GenBank Accession Number:

BC114569

GeneID (NCBI):

54972

UNIPROT ID:

Q24JP5

Full Name:

transmembrane protein 132A

Calculated MW:

1023 aa, 110 kDa

Observed MW:

100-140 kDa

Purification Method:

Protein A purification

CloneNo.:

1H7G8

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : LNCaP cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68589-1-Ig (TMEM132A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.