

# TMEM132A Monoclonal antibody

Catalog Number: 68589-1-Ig

## Basic Information

Catalog Number:

68589-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18149

GenBank Accession Number:

BC114569

GeneID (NCBI):

54972

UNIPROT ID:

Q24JP5

Full Name:

transmembrane protein 132A

Calculated MW:

1023 aa, 110 kDa

Observed MW:

100-140 kDa

Purification Method:

Protein A purification

CloneNo.:

1H7G8

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : LNCaP cells,

## Background Information

### Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68589-1-Ig (TMEM132A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.