

For Research Use Only

# AEG-1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68592-4-PBS



## Basic Information

Catalog Number:

68592-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5033

GenBank Accession Number:

BC045642

GeneID (NCBI):

92140

UNIPROT ID:

Q86UE4

Full Name:

metadherin

Calculated MW:

582 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B11F8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50732-2, AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68592-2-PBS. Detection antibody: 68592-4-PBS. Standard:Ag5033. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50732-4, AEG-1/MTDH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68592-6-PBS. Detection antibody: 68592-4-PBS. Standard:Ag5033. Range: 0.781-100 ng/mL