

For Research Use Only

SLC1A4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68601-2-PBS

Basic Information

Catalog Number:

68601-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17328

GenBank Accession Number:

BC026216

GeneID (NCBI):

6509

UNIPROT ID:

P43007

Full Name:

solute carrier family 1
(glutamate/neutral amino acid
transporter), member 4

Calculated MW:

532 aa, 56 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D1A4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50199-1, SLC1A4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68601-2-PBS. Detection antibody: 68601-3-PBS. Standard:Ag17328. Range: 0.098-12.5 ng/mL