

For Research Use Only

TCF4 Monoclonal antibody

Catalog Number: 68607-1-Ig



Basic Information

Catalog Number:

68607-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17897

GenBank Accession Number:

BC125084

GeneID (NCBI):

6925

UNIPROT ID:

P15884

Full Name:

transcription factor 4

Calculated MW:

671 aa, 72 kDa

Observed MW:

72 kDa; 68 kDa

Purification Method:

Protein G purification

CloneNo.:

1C9A1

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, rat

Positive Controls:

WB : U2OS cells, HCT 116 cells, Caco-2 cells, A549 cells, PC-12 cells, SH-SY5Y cells, K-562 cells, Daudi cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68607-1-Ig (TCF4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.