

For Research Use Only

TCF4 Monoclonal antibody, PBS Only

Catalog Number: 68607-1-PBS



Basic Information

Catalog Number:

68607-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17897

GenBank Accession Number:

BC125084

GeneID (NCBI):

6925

UNIPROT ID:

P15884

Full Name:

transcription factor 4

Calculated MW:

671 aa, 72 kDa

Observed MW:

72 kDa; 68 kDa

Purification Method:

Protein G purification

CloneNo.:

1C9A1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68607-1-Ig (TCF4 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68607-1-PBS in a different storage buffer formulation.