

For Research Use Only

ACSS1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68608-2-PBS



Basic Information

Catalog Number:

68608-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10853

GenBank Accession Number:

BC039261

GeneID (NCBI):

84532

UNIPROT ID:

Q9NUB1

Full Name:

acyl-CoA synthetase short-chain family member 1

Calculated MW:

689 aa, 75 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A9E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50764-1, ACS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68608-2-PBS. Detection antibody: 68608-3-PBS. Standard:Ag10853. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50764-2, ACS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68608-2-PBS. Detection antibody: 68608-4-PBS. Standard:Ag10853. Range: 1.563-100 ng/mL