

For Research Use Only

# TXNDC5 Monoclonal antibody, PBS Only

Catalog Number: 68609-1-PBS



## Basic Information

Catalog Number:

68609-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13965

GenBank Accession Number:

BC001199

GeneID (NCBI):

81567

UNIPROT ID:

Q8NBS9

Full Name:

thioredoxin domain containing 5  
(endoplasmic reticulum)

Calculated MW:

432 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

1A12C3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68609-1-Ig (TXNDC5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68609-1-PBS in a different storage buffer formulation.