

For Research Use Only

UNR; CSDE1 Monoclonal antibody, PBS Only



Catalog Number: 68629-1-PBS

Basic Information

Catalog Number:

68629-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4349

GenBank Accession Number:

BC032446

GeneID (NCBI):

7812

UNIPROT ID:

O75534

Full Name:

cold shock domain containing E1,
RNA-binding

Calculated MW:

798 aa, 89 kDa

Observed MW:

85-90 kDa

Purification Method:

Protein G purification

CloneNo.:

1D9C5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68629-1-Ig (UNR; CSDE1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68629-1-PBS in a different storage buffer formulation.