

For Research Use Only

# TFEB Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68632-2-PBS

## Basic Information

Catalog Number:

68632-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24850

GenBank Accession Number:

BC032448

GeneID (NCBI):

7942

UNIPROT ID:

P19484

Full Name:

transcription factor EB

Calculated MW:

476 aa, 53 kDa

Purification Method:

Protein G purification

CloneNo.:

1D1B9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array in cell lysate using MP50099-1, TFEB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68632-1-PBS. Detection antibody: 68632-2-PBS. Cell lysate: Raji, HepG2, HEK293.



Cytometric bead array standard curve of MP50099-1, TFEB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68632-1-PBS. Detection antibody: 68632-2-PBS. Standard: Ag24850. Range: 0.195-50 ng/mL.