

For Research Use Only

TFEB Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68632-4-PBS

Basic Information

Catalog Number: 68632-4-PBS	GenBank Accession Number: BC032448	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7942	CloneNo.: 1B9C6
Source: Mouse	UNIPROT ID: P19484	
Isotype: IgG1	Full Name: transcription factor EB	
Immunogen Catalog Number: AG24850	Calculated MW: 476 aa, 53 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50099-2, TFEB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68632-3-PBS. Detection antibody: 68632-4-PBS. Standard:Ag24850. Range: 0.781-100 ng/mL.