

For Research Use Only

ST6GAL1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68633-4-PBS



Basic Information

Catalog Number:

68633-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG6181

GenBank Accession Number:

BC040009

GeneID (NCBI):

6480

UNIPROT ID:

P15907

Full Name:

ST6 beta-galactosamide alpha-2,6-sialyltransferase 1

Calculated MW:

47 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2C4H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51249-2, ST6GAL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68633-4-PBS. Detection antibody: 68633-3-PBS. Standard:Ag6181. Range: 0.098-100 ng/mL.