

For Research Use Only

NFE2L1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68634-2-PBS



Basic Information

Catalog Number: 68634-2-PBS	GenBank Accession Number: BC010623	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4779	CloneNo.: 1C4E3
Source: Mouse	UNIPROT ID: Q14494	
Isotype: IgG1	Full Name: nuclear factor (erythroid-derived 2)-like 1	
Immunogen Catalog Number: AG32992	Calculated MW: 772 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50921-1, NFE2L1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68634-2-PBS. Detection antibody: 68634-3-PBS. Standard:Ag32992. Range: 0.195-50 ng/mL