

For Research Use Only

DPF3 Monoclonal antibody

Catalog Number: 68646-1-Ig



Basic Information

Catalog Number:

68646-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG33877

GenBank Accession Number:

BC060801

GeneID (NCBI):

8110

UNIPROT ID:

Q92784

Full Name:

D4, zinc and double PHD fingers, family 3

Calculated MW:

357 aa, 40 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

1F3H4

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, rabbit, rat, mouse, pig

Positive Controls:

WB : pig brain tissue, Jurkat cells, rabbit brain tissue, rat brain tissue, mouse brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68646-1-Ig (DPF3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.