

For Research Use Only

CEACAM1/CD66a Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68653-4-PBS



Basic Information

Catalog Number:

68653-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC014473

GeneID (NCBI):

634

UNIPROT ID:

P13688

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 1 (biliary
glycoprotein)

Calculated MW:

526 aa, 58 kDa

Purification Method:

Protein G purification

CloneNo.:

1D2A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50132-1, CEACAM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68653-3-PBS. Detection antibody: 68653-4-PBS. Standard:Eg0314. Range: 3.125-50 ng/mL.