

For Research Use Only

IL-6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68659-13-PBS



Basic Information

Catalog Number:

68659-13-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1019

GenBank Accession Number:

GeneID (NCBI):

3569

ENSEMBL Gene ID:

ENSG00000136244

Full Name:

interleukin 6 (interferon, beta 2)

Purification Method:

Protein G purification

CloneNo.:

6D1E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50032-3, IL-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68659-13-PBS. Detection antibody: 68659-10-PBS. Standard:HZ-1019 . Range:3.125-50 ng/mL