

For Research Use Only

HGF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68660-1-PBS



Basic Information

Catalog Number:

68660-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

HZ-1084

GenBank Accession Number:

GeneID (NCBI):

3082

Full Name:

hepatocyte growth factor
(hepapoietin A; scatter factor)

Purification Method:

Protein A purification

CloneNo.:

2G2C3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50034-1, HGF Monoclonal Matched Antibody Pair - PBS only. 68660-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1084. 68660-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



Serum of sixteen individual healthy human donors was measured. The human HGF concentration of detected samples was determined to be 2.38 ng/mL with a range of ND - 4.75 ng/mL.