

For Research Use Only

# IL-17F Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68666-1-PBS



## Basic Information

Catalog Number:

68666-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1116

GenBank Accession Number:

GeneID (NCBI):

112744

Full Name:

interleukin 17F

Purification Method:

Protein G purification

CloneNo.:

2C7D11

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50006-1, IL-17F Monoclonal Matched Antibody Pair, PBS Only. 68666-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1116. 68666-2-PBS was HRP conjugated as the detection antibody Range: 62.5-4000 pg/mL