

For Research Use Only

IL-12 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68670-1-PBS

Basic Information

Catalog Number:

68670-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

HZ-1256

GenBank Accession Number:

GeneID (NCBI):

3592

Full Name:

interleukin 12A (natural killer cell
stimulatory factor 1, cytotoxic
lymphocyte maturation factor 1, p35)

Purification Method:

Protein A purification

CloneNo.:

2G1B3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50048-1, IL-12 p40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68670-1-PBS. Detection antibody: 68670-2-PBS. Standard: HZ-1256. Range: 31.25-2000 pg/mL