

For Research Use Only

# PCSK9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68673-2-PBS



## Basic Information

Catalog Number:

68673-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

NM\_174936.4

GeneID (NCBI):

255738

UNIPROT ID:

Q8NBP7

Full Name:

proprotein convertase  
subtilisin/kexin type 9

Calculated MW:

74 kDa

Purification Method:

Protein A purification

CloneNo.:

3H8H9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50055-1, PCSK9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68673-1-PBS. Detection antibody: 68673-2-PBS. Standard: Eg0199. Range: 31.25-2000 pg/mL