

For Research Use Only

FCGR3B/CD16b Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68677-3-PBS



Basic Information

Catalog Number:

68677-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_000570

GeneID (NCBI):

2215

UNIPROT ID:

O75015

Full Name:

Fc fragment of IgG, low affinity IIIb,
receptor (CD16b)

Calculated MW:

26 kDa

Purification Method:

Protein G purification

CloneNo.:

2D9G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50059-2, FCGR3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68677-2-PBS. Detection antibody: 68677-3-PBS. Standard:Eg0308. Range: 0.195-100 ng/mL.