

For Research Use Only

# CD126/IL-6R alpha Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68678-1-PBS



## Basic Information

Catalog Number:

68678-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC132684

GeneID (NCBI):

3570

UNIPROT ID:

P08887

Full Name:

interleukin 6 receptor

Calculated MW:

468 aa, 52 kDa

Purification Method:

Protein G purification

CloneNo.:

2B2A3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50060-1, CD126/IL-6R alpha Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68678-1-PBS. Detection antibody: 68678-2-PBS. Standard: Eg0180. Range: 31.25-2000 pg/mL.