

For Research Use Only

Adiponectin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68679-2-PBS



Basic Information

Catalog Number:

68679-2-PBS

Concentration:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC096308

GeneID (NCBI):

9370

UNIPROT ID:

Q15848

Full Name:

adiponectin, C1Q and collagen domain containing

Calculated MW:

244 aa, 26 kDa

Purification Method:

Protein G purification

CloneNo.:

1F12G8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50050-1, Adiponectin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68679-1-PBS. Detection antibody: 68679-2-PBS. Standard: Eg0327. Range: 62.5-4000 pg/mL.