

For Research Use Only

TNFSF15 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68681-2-PBS



Basic Information

Catalog Number: 68681-2-PBS	GenBank Accession Number: BC074941	Purification Method: Protein A purification
Concentration: 1mg/ml	GeneID (NCBI): 9966	CloneNo.: 1B9F3
Source: Mouse	UNIPROT ID: O95150	
Isotype: IgG2a	Full Name: tumor necrosis factor (ligand) superfamily, member 15	
Immunogen Catalog Number: AG30004	Calculated MW: 28 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50070-1, TNFSF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68681-1-PBS. Detection antibody: 68681-2-PBS. Standard: Ag30004. Range: 78.1-5000 pg/mL.



Serum of sixteen individual healthy human donors was measured. The TNFSF15 concentration of detected samples was determined to be 3,657.9 pg/mL with a range of 2,074.7-4,791.6 pg/mL.